

## IPS e.max® Ceram Impulse Inter Incisal

Inter Incisal is a special incisal material, which is used to increase the brightness value in the incisal third or if the brightness value of the incisal area is high and space is limited.



### Delivery forms

IPS e.max® Ceram Impulse Inter Incisal	1 x 20 g	
White-blue		597006

## IPS e.max® Ceram Impulse Cervical Transpa

With the shaded Cervical Transpa materials, a lifelike in-depth effect is achieved in the cervical third. Compared to the conventional Transpa materials, these special materials demonstrate slightly higher fluorescence.



### Delivery forms

IPS e.max® Ceram Impulse Cervical Transpa	1 x 20 g each	
Yellow		597008
Orange-pink		597009
Kkhaki		597010
Orange		597011

## IPS e.max® Ceram Gingiva

A total of 12 Gingiva ceramic materials for a lifelike gingiva, particularly for the fabrication of implant-supported restorations. Five lifelike gingival materials that are used in combination with the four Intensive Gingiva and three Essence powders, depending on the desires of the patient. An IPS e.max® Ceram ZirLiner Gingiva is available for use in conjunction with zirconium oxide.



### Delivery forms

IPS e.max® Ceram ZirLiner Gingiva	1 x 5 g	605475
IPS e.max® Ceram Essence	1 x 5 g each	
20 rose		605477
21 berry		605478
22 aubergine		605479
IPS e.max® Ceram Gingiva	1 x 20 g each	
1		605481
2		605482
3		605483
4		605484
5		605485
BG34		681357
IPS e.max® Ceram Intensive Gingiva	1 x 20 g each	
1		605487
2		605488
3		605489
4		605490
5		681358