

IPS e.max® Ceram ZirLiner Build-Up Liquid

The IPS e.max® Ceram ZirLiner Build-Up Liquid allround is exclusively used to mix the IPS e.max Ceram ZirLiner powders. The ZirLiner Liquid provides the materials with a creamy and ductile consistency and thus enables easy and efficient application of the ZirLiner. If another consistency is desired, the liquid may also be mixed with the other IPS e.max Ceram Liquids in any mixing ratio.



Delivery forms

IPS e.max® Ceram ZirLiner Build-Up Liquid	1 x 60 ml	
Allround		597050

IPS e.max® Ceram Margin Build-Up Liquids

The IPS e.max® Ceram Margin Build-Up Liquids are available in 2 versions.

- With the allround liquid, a consistency suitable for conventional margin processing and drying is achieved.
- Powders mixed with the carving liquid demonstrate a shorter processing time. After drying, the powder sets. The applied shoulder may then be given the desired shape, i.e. “carved”, using a suitable instrument.



Delivery forms

IPS e.max® Ceram Margin Build-Up Liquid	1 x 60 ml each	
Allround		597052
Carving		597053