

IPS® Build-Up Liquid

The IPS® Build-Up Liquids are suitable for both IPS e.max® Ceram and IPS Style® Ceram. They are available in two versions:

- With the allround liquid, a consistency suitable for the conventional processing and drying is achieved.
- The soft liquid demonstrates shorter processing times and faster drying.



Delivery forms

IPS® Build-Up Liquid	1 x 60 ml
Allround	673375
Soft	673376
IPS® Build-Up Liquid	1 x 250 ml
Allround	673377
Soft	673378

IPS e.max® Ceram Glaze and Stain Liquid

The IPS e.max® Ceram Glaze and Stain Liquids are available in 2 versions:

- With the allround liquid, a consistency suitable for conventional processing and drying is achieved. Powders mixed with the all-round liquid demonstrate a shorter processing time (approximately 1/2 day).
- The longlife liquid, however, enables a slightly pasty consistency and the powders mixed with this liquid can be processed over a longer period of time (approximately 1 week).



Delivery forms

IPS e.max® Ceram Glaze and Stain Liquid	1 x 15 ml each
Allround	597058
Longlife	597059