

## IPS® Multi Wax Pattern Form A/Form B

The IPS® Multi Wax Patterns are precision wax patterns suitable for spruing the wax-up to the ring base. The special shape of the wax patterns decisively influences the flow behaviour of the IPS e.max Press Multi material during the press procedure. The patterns ensure that the progression of shade (dentin/incisal) is pressed into the correct position in the restoration. Form B is used to sprue very delicate restorations (e.g. mandibular anterior teeth) to the ring base. Form A is used for all the other types of restorations.



### Delivery forms

<b>IPS® Multi Wax Pattern Form A</b>	<b>1 x 80</b>	<b>638163</b>
<b>IPS® Multi Wax Pattern Form B</b>	<b>1 x 80</b>	<b>638164</b>

## IPS® Multi Investment Ring System 200 g

The IPS® Multi Investment Ring Base 200 g is used in the processing of IPS e.max® Press Multi ingots. The ring base features openings which can accommodate up to four wax-ups. The IPS Multi Investment Rings Base 200 g is offered together with the IPS Ring Gauge and the IPS Silicone Ring as the IPS Multi Investment Ring System Kit 200 g.



### Delivery forms

<b>IPS® Multi Investment Ring System 200 g</b>	<b>3x</b>	<b>638606</b>
--	-----------	---------------