

IPS® Alox Plunger Separator

The IPS® Alox Plunger Separator prevents the ingot from sticking to the Alox plunger during pressing and cooling in the investment ring by forming a separation layer. A very small amount of separator is sufficient.



Delivery forms

IPS® Alox Plunger Separator	1 x 200 mg	597063
-----------------------------	------------	--------

IPS® Ingot Tong

The IPS® Ingot Tong are used to simultaneously place several press ingots in the 300-g investment ring.



Delivery forms

IPS® Ingot Tong	1x	614498
-----------------	----	--------

IPS e.max® Press Invex Liquid

The IPS e.max® Press Invex (Investment-Ex) Liquid is used to soften the surface reaction layer formed on pressed IPS e.max Press and ZirPress objects during the pressing procedure.



Delivery forms

IPS e.max® Press Invex Liquid	1 x 1 l	597064
IPS e.max® Press Invex jar (empty)	1x	628336
IPS e.max® Press Invex Liquid	1 x 50 ml	630539