

IPS e.max® Press / IPS e.max® CAD Material Shade Guides

The shade guides for press and CAD materials are supporting tools for determining the shade. To prevent cross-contamination in the oral cavity, the individual tabs and the entire shade guide are autoclavable. The shade tabs show the final shade, i.e. after pressing or crystallization. They thus help to check the shade, and they can also be used to check the quality of the crystallization process of IPS e.max® CAD restorations.



Delivery forms

IPS e.max® Press HO/MO Shade Guide	1x	597065
IPS e.max® CAD MO Shade Guide	1x	597067
IPS e.max® Press / IPS e.max CAD LT Shade Guide	1x	605523
IPS e.max® Press / IPS e.max CAD HT Shade Guide	1x	613998
IPS e.max® Press / IPS e.max CAD MT Shade Guide	1x	680311

IPS e.max® Ceram Material Shade Guide

For IPS e.max® Ceram, the Material Shade Guides are an integral part of the product concept. To prevent cross-contamination in the oral cavity, the individual tabs and the entire shade guide are autoclavable. For reasons of light refraction, the dentin shade tabs have been given an anatomical shape and surface structure. All the other shade tabs demonstrate a slight wedge shape, in order to better illustrate the translucency of the individual materials. To facilitate the distinction from the existing material shade guides (e.g. IPS d.SIGN®, IPS InLine®), the rods and holders of the IPS e.max Ceram material shade guides are light-grey.



Delivery forms

IPS e.max® Ceram Material Shade Guide	1 x each	
IPS e.max Ceram ZirLiner		597069
IPS e.max Ceram Margin A–D		597071
IPS e.max Ceram Deep Dentin A–D		597073
IPS e.max Ceram Power Dentin/Incisal		681621
IPS e.max Ceram Dentin A–D		597276
IPS e.max Ceram Incisal/Transpa		597074
IPS e.max Ceram Impulse		597075
IPS e.max Ceram Essence		597076
IPS e.max Ceram Bleach BL		602463
IPS e.max Ceram Gingiva		602464
IPS e.max Ceram Selection		685944

IPS e.max® Packaging

The IPS e.max® Material Box and the inserts are the ideal storage place in the laboratory for the various IPS e.max materials.



Delivery forms

IPS e.max® Packaging		
IPS e.max Material Box small (55 mm)	1x	597483
IPS e.max Material Box medium (80 mm)	1x	597484
IPS e.max Ceram Insert (for 22 refills)	1x	597485
IPS e.max Ceram Shade Insert (for 16 refills + 2 liquids)	1x	635843
IPS e.max Ceram Essence Insert (for 21 refills + 2 liquids)	1x	635844
IPS e.max Press/ZirPress Insert (for 24 refills)	1x	607668
IPS e.max CAD/ZirCAD Insert (for 6 refills)	1x	614248
IPS e.max Press Multi Ingot Insert (for 24 ingots)	1x	672835