

IPS Empress® Esthetic Veneer Wash Pastes

With the six standard Wash pastes, the pressed/reduced restorations are individually characterized. The ideal gradation of shade nuances allows the fabrication of highly esthetic restorations. Moreover, the five Wash paste modifiers (copper, orange, sky blue, basic yellow, basic red) permit the individual shading of the six standard Wash pastes.



Delivery forms

| IPS Empress® Esthetic Veneer Wash Pastes | 1 x 1 g each | |
|--|--------------|--------|
| Neutral | | 579406 |
| Low value | | 579407 |
| High value | | 579408 |
| MM light-salmon | | 579409 |
| MM yellow-orange | | 579410 |
| MM reddish-orange | | 579411 |
| Modifier copper | | 579412 |
| Modifier orange | | 579413 |
| Modifier sky blue | | 579414 |
| Modifier basic yellow | | 579415 |
| Modifier basic red | | 579416 |

IPS Empress® Esthetic Veneer Incisal Opal

These opalescent materials are available in three different gradations (low, medium, and high translucency).



Delivery forms

| IPS Empress® Esthetic Veneer Incisal Opal | 1 x 20 g each | |
|---|---------------|--------|
| Low translucent | | 579388 |
| Medium translucent | | 579389 |
| High translucent | | 579390 |

IPS Empress® Esthetic Veneer Incisal

The Incisal materials are available in six shade nuances (white, grey, orange, edge, yellow, orange-pink) and are used to apply additional characterizations in the incisal area.



Delivery forms

| IPS Empress® Esthetic Veneer Incisal | 1 x 20 g each | |
|--------------------------------------|---------------|--------|
| White | | 579391 |
| Grey | | 579392 |
| Orange | | 579393 |
| Edge | | 579394 |
| Yellow | | 579397 |
| Orange-pink | | 579398 |

IPS Empress® Esthetic Veneer Chroma Modifier

Supports the shade enhancement of the pressed/reduced restoration and is available in the shades 110/A1 and 210/A3.



Delivery forms

| IPS Empress® Esthetic Veneer Chroma Modifier | 1 x 20 g each | |
|--|---------------|--------|
| 110/A1 | | 579395 |
| 210/A3 | | 579396 |