

IPS Empress® Esthetic Veneer Brightener

To increase the brightness value



Delivery forms

IPS Empress® Esthetic Veneer Brightener	1 x 20 g	579399
-----------------------------------------	----------	--------

IPS Empress® Esthetic Veneer Transparent

The shaded Transparent materials are available in two nuances (blue and neutral).



Delivery forms

IPS Empress® Esthetic Veneer Transpa	1 x 20 g each	
Blue		579400
Neutral		579401

IPS Empress® Add-On 770 °C / 1418 °F

Low-fusing add-on material for subsequent adjustments, such as proximal and occlusal contacts for:

- IPS Empress® Esthetic restorations
- IPS Empress layering technique restorations



Delivery forms

IPS Empress® Add-On 770 °C / 1418 °F	1 x 20 g	564080
--------------------------------------	----------	--------

IPS Empress® Esthetic Veneer Build-Up Liquid

This liquid is used to mix the IPS Empress® Esthetic Veneer layering materials and IPS Empress Add-On 770°C / 1418°F.



Delivery forms

IPS Empress® Esthetic Veneer Build-Up Liquid	1 x 60 ml	579405
----------------------------------------------	-----------	--------