

## IPS Empress® Investment Ring System

The Investment Ring System is used to sprue the contoured restorations. The 100-g Investment Ring System is suitable for single restorations, while the 200-g Investment Ring System is used for several single restorations. Given the different diameters of the pressed ingots and Alox plungers, the IPS Empress® Investment Ring System is only suitable for the materials of the IPS Empress System.



### Delivery forms

IPS Empress® Investment Ring System		
Investment Base	3x	
Investment Gauge	3x	
100 g		<b>554034</b>
200 g		<b>554035</b>

## IPS® Silicone Ring

- IPS® Silicone Ring 100 g for the IPS Empress® and IPS e.max® Investment Ring System 100 g
- IPS Silicone Ring 200 g for the IPS Empress and IPS e.max Investment Ring System 200 g



### Delivery forms

IPS® Silicone Ring	1 x each	
100 g		<b>590113AN</b>
200 g		<b>590114AN</b>

## IPS Empress® Alox Plunger

With the IPS Empress® Alox Plunger, the ceramic ingot in its softened state is pressed into the hollow space in the investment ring. It features a slightly smaller diameter and is somewhat longer than the IPS e.max® Alox Plunger and is, therefore, exclusively suitable for the IPS Empress System.



### Delivery forms

IPS Empress® Alox Plunger	1x	<b>554070</b>
---------------------------	----	---------------