

IPS® PressVEST Premium

IPS® PressVEST Premium is a universal, phosphate-bonded investment material for the following Ivoclar Vivadent press ceramics:

- IPS e.max® Press
- IPS e.max ZirPress
- IPS InLine® PoM (Press-on-Metal ceramic)
- IPS Empress Esthetic

IPS PressVest Premium can be used in the popular speed heating procedure as well as with conventional heating methods. The high stability during the preheating phase allows materials with different expansion values to be processed. IPS PressVest Premium is therefore ideally suitable for milled and printed objects made from wax or acrylic. Due to the material's optimum properties, IPS PressVest Premium has only a slight reaction layer, which is very easy to remove, and allows for a silky smooth, homogeneous surface.

IPS PressVEST Premium **cannot** be used for metal casting, due to the high temperatures.

Delivery forms

IPS® PressVEST Premium

Powder	1 x 2.5 kg	685585AN
	1 x 5 kg	685586AN
Liquid	1 x 500 ml	685587
	1 x 1 l	685588

