

IPS® Investment Ring Tongs

The IPS® Investment Ring Tongs have been designed to ensure safe working procedures. The tongs are used to place the press ceramic ingots and Alox plunger in the investment ring prior to pressing.



Delivery forms

IPS® Investment Ring Tongs	1x	554073
----------------------------	----	--------

Tweezers



Delivery forms

Tweezers	1x	554069
----------	----	--------

IPS® Spatula



Delivery forms

IPS® Spatula	1x	635529
--------------	----	--------

Modelling Brushes

These red marten modelling brushes are used to layer the ceramic materials.

Sizes: 00, 0, 1, 2, 4, 6, 10, G2, G4



Delivery forms

Modelling Brushes Assortment		530381
-------------------------------------	--	---------------

Modelling Brushes (sizes: 00, 0, 1, 2, 4, 6, 10, G2, G4)	9x	
--	----	--

Modelling Brushes Refill	1x each	
---------------------------------	----------------	--

00		530382
0		530383
1		530384
2		530385
4		530386
6		530387
10		530388
G2		530389
G4		530390