

# IPS® Natural Die Material Shade Guide

Determination of the preparation shade in the fabrication of all-ceramic restorations

The IPS Natural Die Material shade guide enables optimum selection of the shade after preparation. The shade is selected on the basis of the most severely discoloured part of the prepared tooth.

The IPS Natural Die Material allows you to fabricate a customized model die with which you can then accurately select the shade and brightness of all-ceramic restorations.

## Indication:

- Determination of the preparation shade after tooth preparation

## Advantages:

- Easy and accurate reproduction of shades due to having determined the preparation shade
- Fabrication of a customized model die
- Autoclavable

## Shades:


- 9 (ND 1, ND 2, ND 3, ND 4, ND 5, ND 6, ND 7, ND 8, ND 9)



## Delivery forms

<b>IPS® Natural Die Material Kit</b>		<b>597078</b>
IPS Natural Die Material, all 9 shades	9 x 8 g	
IPS Natural Die Material Separator	1 x 20 g	
IPS Natural Die Material Shade Guide	1x	
Various accessories		
<b>IPS Natural Die Material Kit Promo</b>		
IPS Natural Die Material Kit (597078)	1x	
IPS Natural Die Material Shade Guide	4x	
<b>Refills IPS Natural Die Material</b>	<b>8 g each</b>	
<b>ND 1</b>		<b>597080</b>
<b>ND 2</b>		<b>597081</b>
<b>ND 3</b>		<b>597082</b>
<b>ND 4</b>		<b>597083</b>
<b>ND 5</b>		<b>597084</b>
<b>ND 6</b>		<b>597085</b>
<b>ND 7</b>		<b>597086</b>
<b>ND 8</b>		<b>597087</b>
<b>ND 9</b>		<b>597088</b>
<b>IPS Natural Die Material Separator</b>	<b>1 x 20 g</b>	<b>597089</b>
<b>Die pins</b>	<b>10x</b>	<b>554028</b>
<b>IPS Condenser</b>	<b>10x</b>	<b>554030</b>
<b>IPS Natural Die Material Shade Guide (autoclavable)</b>		<b>597090</b>

## Influence of the preparation shade



Restoration: veneer (0.5 mm, IPS e.max® CAD HT B1)  
Luting material: Variolink® Esthetic, neutral  
Preparation shade: IPS® Natural Die Material, ND 1 – ND 9