

IPS® Object Fix Flow / IPS® Object Fix Putty

IPS® Object Fix auxiliary firing pastes with different viscosity to support all-ceramic restorations during firing and crystallization. The pastes are used to stabilize and secure the restorations on the silicon nitride support pins or directly on the silicon nitride firing tray. IPS Object Fix Putty and IPS Object Fix Flow are easy to apply and remove after the firing procedure.



Delivery forms

IPS® Object Fix

Flow	1 x 10 ml	605525AN
Putty	1 x 10 g	605526 AN

All-Ceramic Preparation Guide

The three-dimensional preparation guide from Ivoclar Vivadent visualizes the proper preparation of teeth for all-ceramic restorations. The model furthermore shows the minimum preparation requirements and the indications of IPS Empress® and IPS e.max® restorations. The model comes with an illustrated booklet which features the preparation recommendations. The booklet is also available separately.



Delivery forms

All-Ceramic Preparation Guide Edition II	630548
All-Ceramic Preparation Guidelines (D/E)	631797