

IPS Ivocolor® Starter Kit

The IPS Ivocolor® Starter Kit is the universal stains and glaze range for all IPS ceramic materials. The IPS Ivocolor Starter Kit features the most popular shade pastes, essence powders, a glaze paste, a glaze powder, two mixing liquids and two material shade guides. The IPS Ivocolor Starter Kit is supplied in the material box and can be individually supplemented with other products.



Delivery forms

| IPS Ivocolor® Starter Kit | 667649 |
|--|----------------|
| IPS Ivocolor Essence, shades: E 01 white, E 02 creme, E 10 mahogany, E 17 anthracite, E 21 basic red | 1 x 1.8 g each |
| IPS Ivocolor Shade Dentin, shades: SD1, SD2, SD4 | 1 x 3 g each |
| IPS Ivocolor Shade Incisal, shades: SI1, SI3 | 1 x 3 g each |
| IPS Ivocolor Glaze Powder | 1x 1.8 g |
| IPS Ivocolor Glaze Paste | 1 x 3 g |
| IPS Ivocolor Mixing Liquid (allround) | 1 x 15 ml |
| IPS Ivocolor Essence Fluid | 1 x 15 ml |
| IPS Ivocolor Essence Material Shade Guide | 1x |
| IPS Ivocolor Shade Material Shade Guide | 1x |
| Modelling brushes | 1x |

IPS Ivocolor® Essence Kit

The IPS Ivocolor® Essence Kit contains all IPS Ivocolor Essence shade powders, glaze powders and the corresponding mixing liquids. The IPS Ivocolor Essence Kit comes with the Essence Material Shade Guide.



Delivery forms

| IPS Ivocolor® Essence Kit | 667650 |
|--|----------------|
| IPS Ivocolor Essence, shades: E 01 white, E 02 creme, E 03 lemon, E 10 mahogany, E 17 anthracite, E 21 basic red, E 04 sunset, E 05 copper, E 06 hazel, E 07 olive, E 08 khaki, E 09 terracotta, E 10 mahogany, E 11 cappuccino, E 12 espresso, E 13 terra, E 14 profundo, E 15 ocean, E 16 sapphire, E 17 anthracite, E 18 black, E 19 rose, E 20 coral, E 21 basic red, E 22 basic yellow, E 23 basic blue | 1 x 1.8 g each |
| IPS Ivocolor Glaze Powder | 1 x 5 g |
| IPS Ivocolor Glaze Powder FLUO | 1 x 5 g |
| IPS Ivocolor Mixing Liquid (allround) | 1 x 15 ml |
| IPS Ivocolor Essence Fluid | 1 x 15 ml |
| IPS Ivocolor Essence Material Shade Guide | 1x |

IPS Ivocolor® Shade Kit

The IPS Ivocolor® Shade Kit contains nine Shade Dentin materials, three Ivocolor Shade Incisal materials, glaze pastes and mixing liquids. The kit is supplemented by the IPS Ivocolor Shade Material Shade Guide.



Delivery forms

| IPS Ivocolor® Shade Kit | 667651 |
|--|--------------|
| IPS Ivocolor Shade Dentin, shades: SD0, SD1, SD2, SD3, SD4, SD5, SD6, SD7, SD8 | 1 x 3 g each |
| IPS Ivocolor Shade Incisal, shades: SI1, SI2, SI3 | 1 x 3 g each |
| IPS Ivocolor Glaze Paste | 1 x 3 g |
| IPS Ivocolor Glaze Paste FLUO | 1 x 3 g |
| IPS Ivocolor Mixing Liquid (allround) | 1 x 15 ml |
| IPS Ivocolor Mixing Liquid (longlife) | 1 x 15 ml |
| IPS Ivocolor Shade Material Shade Guide | |