

IPS Ivocolor® Glaze

- **IPS Ivocolor® Glaze Paste** is a new, thixotropic glaze paste for a high degree of gloss in only one firing cycle. The glaze is available with or without fluorescence and is compatible with all IPS ceramics.
- **IPS Ivocolor Glaze Powder** impresses users with its excellent mixing and handling behaviour and is ideally suited for the glazing of layering ceramics.



Delivery forms

IPS Ivocolor® Glaze Powder	1 x 1.8 g	667686
	1 x 5 g	667688
IPS Ivocolor® Glaze Powder FLUO	1 x 1.8 g	667687
	1 x 5 g	667689
IPS Ivocolor® Glaze Paste	1 x 3 g	667690
	1 x 9 g	667692
IPS Ivocolor® Glaze Paste FLUO	1 x 3 g	667691
	1 x 9 g	667693

IPS Ivocolor® Mixing Liquids

IPS Ivocolor® Mixing Liquids are suitable for the mixing and diluting of all IPS Ivocolor powders and pastes.

- **IPS Ivocolor Mixing Liquid allround** has a low viscosity. Mixed powders remain applicable for a shorter period of time (approx. ½ day).
- **IPS Ivocolor Mixing Liquid longlife** produces a higher viscosity and remains applicable for a longer period of time (approx. 1 week).



Delivery forms

IPS Ivocolor® Mixing Liquid allround	1 x 15 ml	667694
IPS Ivocolor® Mixing Liquid longlife	1 x 15 ml	667695

IPS Ivocolor® Essence Fluid

IPS Ivocolor® Essence Fluid is only suitable for the mixing of IPS Ivocolor Essence. Mixing IPS Ivocolor Essence Fluid with the shade powders results in a paste-like consistency.



Delivery forms

IPS Ivocolor® Essence Fluid	1 x 15 ml	667696
-----------------------------	-----------	--------