

## Vacuum Pump Filter Kit

This vacuum pump filter kit protects the vacuum pump from contamination.



### Delivery forms

<b>Vacuum Pump Filter Kit</b>		<b>680188</b>
Vacuum Pump Filter Kit		1x
Vacuum hose 200 mm		1x

## Programat® Firing Tray Kit 2

The Firing Tray Kit 2 was specially developed for the Programat® furnaces. The support pins can be individually mounted to the firing tray if required.



### Delivery forms

<b>Programat® Firing Tray Kit</b>		<b>655946</b>
Firing trays		2x
Support pins A		6x
Support pins B		6x
Fan pins		2x
Ceramic sleeves		2x
Hobbyceram pen		1x

## Metal Pins for Firing Tray

The support pins can be used for the Programat® firing tray and can be mounted to the honeycomb firing tray if required.



### Delivery forms

<b>Metal Pins for Firing Tray, Refill</b>		<b>626698</b>
Metal pins A		4x
Metal pins B		4x
Metal pins C		4x

## Ceramic Sleeves for Metal Pins

The ceramic sleeves can be placed on the metal pins to avoid that metal-free restorations come in contact with the metal pins.



### Delivery forms

<b>Ceramic Sleeves for Metal Pins Refill</b>		<b>666068</b>
Ceramic sleeves		10x