

Firing Tongs

These firing tongs are used to place objects in the firing chamber. High precision and excellent material quality are the special features of these tongs.



Delivery forms

Firing Tongs	1x	530665
---------------------	-----------	---------------

Programat® Tong Holder

The Programat® tong holder can be mounted to the Programat housing and permits convenient placement of the firing tongs.

For Programat P300, P500, P700,
EP 3000, EP 5000



For Programat CS2, CS3, P310, P510,
P710, EP 3010, EP 5010



Delivery forms

Programat® Tong Holder	1x	607538
-------------------------------	-----------	---------------

(for Programat P300, P500, P700, EP 3000, EP 5000)

Programat® Tong Holder xx10	1x	655945
------------------------------------	-----------	---------------

(for Programat CS2, CS3, P310, P510, P710, EP 3010, EP 5010)

USB Download Cable

You can update/upgrade the equipment software via PC / laptop using the USB download cable. In addition, the USB download cable is suitable to control the furnace by means of the PrograBase 3 or PrograBase X10 software from the PC and to protocol firing data.



Delivery forms

USB Download Cable Type A-B, 1.8 m	1x	597183
---	-----------	---------------