

Programat® Sinter Tray

The sinter tray permits easy loading of the Programat® S1/ Programat S1 1600 and is imperative for optimum sinter results. The special geometry of the tray ensures homogeneous heat distribution.



Delivery forms

Programat® Sinter Tray	1x	627987
------------------------	----	--------

Programat® Dosto Tray

With the new Programat® Dosto Tray, the loading capacity of the Programat S1 and the Programat S1 1600 can be increased by approx. 40% per firing cycle (with the exception of speed programs). This makes sintering even more efficient.



Delivery forms

Programat® Dosto Tray	1x	688790
-----------------------	----	--------

Programat® Sinter Tray Fork

The sinter tray fork is used to place the sinter tray inside the sintering chamber. Its shape has been precisely adjusted to fit the Programat sinter tray.



Delivery forms

Programat® Sinter Tray Fork	1x	627988
-----------------------------	----	--------