



### ★ STARMOD

The ideal modeling wax for gnathological waxing techniques. Available in 2 different hardnesses: Hard / 57°C, Regular / 55°C.

- › High stability
- › Precise drawing of the wax
- › Leaves smooth surface

VERSION	HARDNESS	REF
70g grey	hard	521010H
70g ivory	hard	521020H
70g blue	hard	521030H
70g ivory	regular	521020
70g blue	regular	521030



### ★ STARMOD F / STARMOD F PUR

In combined techniques different demands are made of the waxes. In the case of pure milling, nothing should lubricate, the wax must be hard. When scraping on the other hand, the wax should be not too hard. Since it is almost impossible to combine both requirements in one wax, we have adapted our milling waxes to the respective application areas.

- › StarMod F = Milling and Scraping, working temperature 57°C
- › StarMod F Pur = Milling, working temperature 58°C

VERSION	APPLICATION	REF
70g StarMod F green	milling and scraping	521100
70g StarMod F mint green	milling	521130



### ★ DIAGNOSTIC WAXSTAR

A wax which is suitable for an aesthetic wax-up as well as for the pressed ceramic at the same time. By using purely organic ingredients, Diagnostic WaxStar burns absolutely without residue. The Vita basic colors A-B-C-D as well as white ensure a fitting, which is almost the finished work. Diagnostic wax is available in the 70 g tin, or as a complete set in 5 x 9 g.

- › Working temperature: 54 °C



VERSION	REF
70g white	521700
70g dentin A	521710
70g dentin B	521720
70g dentin C	521730
70g dentin D	521740
5x 9g SET with all 5 colour shades	521750



OK  
PRESSABLE  
CERAMIC

VERSION	CROSS-SECTION	REF
70g round	Ø 3mm	522210
150g round	Ø 4mm	522230
150g round	Ø 5mm	522240
70g square	2,5x2,5mm	522110

### ★ STICKYSTAR

The adhesive sticks enable better treatment and prevent the distortion of the casting object when removing from the model.

### ★ SPRUESTAR / SPRUESTAR XL

SprueStar is the optimised wax casting pear, with the result, that less turbulence occurs in the injected melting due to the nozzle shape of the casting duct. A continuous treating of the casting is simplified and improved by flooding with the small waxing casting stick.

- › Simple and metal saving separating
- › Duct diameter only 3mm
- › No blowholes
- › 3 different reservoir sizes

- › SprueStar XL (for NPM)
- › Optimised wax casting pear



OK  
PRESSABLE  
CERAMIC

VERSION	RESERVOIR	REF
300 Sticks with wax casting pear	Ø 6mm	521900
200 Sticks with wax casting	Ø 7mm	521905
155 Casting sticks	Ø 8mm	521910
200 XL Sticks with wax casting pear	Ø 8mm	521880
180 XL Sticks with wax casting pear	Ø 9mm	521890



OK  
PRESSABLE  
CERAMIC

### ★ WIRESTAR HARD

An ideal wire wax for crown, bridge and model casting technique, which can be flexed and deformed without any effort.

- › Dimensionally stable and tension-free

VERSION	Ø	REF
250g green	2,5mm	522010H
250g green	3mm	522020H
250g green	3,5mm	522030H
250g green	4mm	522040H
250g green	5mm	522050H