



✦ O-SPRAY

Occlusion spray for the marking of contact points, areas and impurities. In four easily recognisable colours for optimal contrasts

- › Uniform, fine spray pattern, economical
- › Precise dosage
- › Easy removal with water
- › Reduced dust behavior
- › Food ingredients

VERSION	REF
75ml white	230233
75ml blue	230234
75ml red	230235
75ml green	230236



✦ DIASTAR

All-in-One diamond polished paste for metals, ceramics, zirconia, composites, plastics and acrylics.

- › Very high diamond content
- › Hard paste consistency, thereby more efficient dosing
- › No splashing or smearing - clean and economical
- › Not suitable for intraoral use!

VERSION	REF
1 can 30g	271069



✦ BLACK DIAMOND

Diamond polishing paste for all dental alloys, including NPM and titanium.

- › Does not lubricate
- › Economical in use
- › Perfect polishing results

VERSION	REF
450g	271071



✦ CROWN FORCEPS

Made of flexible, highly heat-resistant, light carbon fiber. Excellent for fixing various works and materials. Auxiliary instrument when separating, finishing and polishing, without hot fingers. Especially perfect for cleaning and evaporation without loss of the workpiece.

VERSION	REF
1 pieces	230240



✦ ACID TWEEZER

Made of teflon-coated steel. The tweezers are used to prevent the metal contact when inserting or removing cast objects in acid / pickling solutions.

VERSION	REF
1 piece	230800



✦ HANDY TORCH LABO

The blowtorch "Handy Torch Labo" is very easy to use. It can be operated hand-held or in standing position at the workstation.

- › Automatic piezo ignition
- › Flame regulation
- › Maximum temperature 1300 °C

VERSION	REF
1 piece	500100

09

MATERIALS FOR THE CERAMIST

Premium brushes for the ceramist and accessories for ceramic firing.

- ✦ DENTAL BRUSH SCHEFTNER KOLINSKY
- ✦ CERAMIC BRUSH SCHEFTNER TORAY
- ✦ ISOKER FLUID / ISOKER SPRAY
- ✦ STARTRAY HONEYCOMB FIRING TRAYS
- ✦ CERAMIC HOLDING PINS SET
- ✦ HOLDING PIN SET
- ✦ PINSSTAR / PINSSTAR UNIVERSAL REVERSIBLE
- ✦ FIBRE TRAY

