

## MPF CROWN HOLDER


NEW

Assistive tool to safely hold crowns in place while working, without chipping or marking.

**Item No. 152-0001**  
MPF Crown Holder 3 - 4mm

**Item No. 152-0002**  
MPF Crown Holder 2 - 1mm

The MPF Crown Holder is an instrument designed to successfully hold crowns in place while working, without chipping or marking them. Special gum material is applied to the desired end (retention knurl or cut out), onto which the crown is positioned. The adhesive gum is guaranteed not to leave sticky, greasy residue on the crown or on fingers.



**Item No. 152-0000** MPF Crown Holder Complete Kit  
Includes 4 - 3mm & 2 - 1mm handles and 1pk. White Gum & 1pk. Black Gum Material

**Item No. 152-0003** MPF White Gum Material  
**Item No. 152-0004** MPF Black Gum Material

## Veneer & Crown Stabilizing Paste

NEW



A stabilizing paste used to safely abrade metal porcelain crowns, all ceramic crowns and veneers, composite crowns and veneers without the prosthesis detaching itself. Burns residue free.

**Item No. 150-0001** Veneer & Crown Stabilizing Paste

## Peg Fix



A new, revolutionary, environmentally friendly refractory stabilizing paste, made from PCW fibers and suitable for firing with zirconia, metal or veneers. Can be completely removed from your prosthesis with a short blast of compressed air, leaving it totally residue free.

**Item No. 150-0002** Peg Fix