

## ZMax

FOR ZIRCONIA

Diamond impregnated discs for zirconia abrasion. Designed to generate minimal heat.

**Item No. 120-0001**  
ZMax Large Disc, 22 x 4.5mm



**Item No. 120-0002**  
ZMax Small Disc, 15 x 3.5mm



**Item No. 120-0003**  
ZMax Disc, 22 x 2mm



**Item No. 120-0004**  
ZMax Thin Disc, 22 x 0.7mm



**Item No. 120-0005**  
ZMax Small Barrel, 4.8 x 13mm



**Item No. 120-0006**  
ZMax Large Barrel, 6.5 x 13mm



**Item No. 120-0007**  
ZMax Cone, 3.5 x 11mm



**Item No. 120-0008**  
ZMax Inverted Cone, 12 x 6mm



**Item No. 120-0009**  
ZMax Inverted Small Cone, 6 x 8mm



## Zirco Cera GSS

FOR ZIRCONIA & CERAMIC

Diamond impregnated discs for both zirconia & ceramic abrasion. Minimal heat & pressure.

**Item No. 121-0001**  
Zirco Cera GSS Large Disc  
22 x 4.5mm



**Item No. 121-0002**  
Zirco Cera GSS Small Disc  
15 x 3.5mm



**Item No. 121-0003**  
Zirco Cera GSS Disc  
22 x 2mm



**Item No. 121-0004**  
Zirco Cera GSS Thin Disc  
22 x 0.7mm



**Item No. 121-0005**  
Zirco Cera GSS Small Barrel  
4.8 x 13mm



**Item No. 121-0006**  
Zirco Cera GSS Large Barrel  
6.5 x 13mm



**Item No. 121-0007**  
Zirco Cera GSS Cone  
3.5 x 11mm



**Item No. 121-0008**  
Zirco Cera GSS Inverted Cone  
12 x 6mm



**Item No. 121-0009**  
Zirco Cera GSS Inverted Small Cone  
6 x 8mm

