

MPF RAINBOW CERAMIC POLISHER

Premium diamond abrasive wheels for perfectly sealed ceramic or zirconia restoration.



Item No. 122-0001 MPF Rainbow Ceramic Polisher
Coarse, Red, RPM 13000-15000

For minimal reduction on zirconia marginal regions and all types of ceramic materials.




Item No. 122-0002 MPF Rainbow Ceramic Polisher
Medium, Blue, RPM 13000-15000

For minimal reduction on zirconia marginal regions and for contouring all types of ceramic materials.



Item No. 122-0003 MPF Rainbow Ceramic Polisher
Fine, Grey, RPM 10000-13000

For the removal of minor scratches on all types of ceramic materials.



Item No. 122-0004 MPF Rainbow Ceramic Polisher
High Shine, Yellow, RPM 8000-10000

For the final high shine procedure on all types of ceramic and zirconia materials.

Polishing Brushes



Item No. 124-0001 Polishing Brush, 100pcs.
Item No. 124-0002 Polishing Brush, 50pcs.
Item No. 124-0003 Polishing Brush, 25pcs.

Made from natural hair. For an exceptional alloy polish.



Item No. 125-0001 Polishing Mop, 100pcs.
Item No. 125-0002 Polishing Mop, 50pcs.
Item No. 125-0003 Polishing Mop, 25pcs.

Ideal thickness & hardness for polishing ceramic surfaces.

Aluminium Oxide Abrasive Disc



Item No. 123-0001 Middle 25 x 0.6 x 1.8mm, 100pcs.
Item No. 123-0002 Large 38 x 0.6 x 1.8mm, 100pcs.
Item No. 123-2003 Small 22 x 0.2-0.25 x 1.8mm, 100pcs.

Highly durable discs for ceramic and alloy restoration, without the risk of chipping or cracking the surface.