

An easy way to make a great impression



VivaScanTM

The compact and intuitive
intraoral scanning solution
for your practice

Simplicity meets performance

VivaScan is a compact and intuitive intraoral scanning solution for your practice. Due to its standalone design the device can be easily and effortlessly integrated into your practice to make it more efficient. VivaScan is the ideal solution for your entry into the world of digital dentistry.


- ✓ High-accuracy 3D scanning for high-precision digital impressions
- ✓ Intuitive software and easy to use
- ✓ Compact in size for a flexible scanning experience
- ✓ Increased patient comfort
- ✓ Direct connection with your laboratory
- ✓ Digital entry to IPS e.max® world- leading all-ceramics

Training & education

For a successful start with your new VivaScan, we offer you web-based, on-demand instructions on how to use the scanning device. We provide additional on-site training and education measures if needed.

Handling made easy

VivaScan enables you to make digital impressions of your patients' teeth quickly, easily and effortlessly. The intraoral scanner is connected to the laptop with only one cable, which contributes to a convenient handling. High speed provides a fast scanning experience, helping you to treat more patients without wasting time or material.



Plug & play function

Allows you to use VivaScan as a portable system after it is connected to a laptop

Two button function

Allows a touch-free operation on the patient

Two tip sizes

Allow you to determine the best fit for every patient

True color scan

Helps your lab to design and fabricate restoration precisely

Integrated HD camera

Makes your work through high-resolution photographs and case documentation easier

Lightweight and ergonomic design

Makes scanning easy for everyone

Ivoclar Cloud

VivaScan's straightforward, user-friendly and intuitive software lets you send scans directly to your laboratory of choice in a streamlined workflow. Our integrated file sharing solution ensures secure transfer of data.

Never been more connected than today

One click

You simply upload your scans to the Ivoclar Cloud with just one click and your technician can download them just as easily.

Open file format

All files can be downloaded from the cloud in an open file format, so your technician is free to work in his or her preferred way.

No subscription fee

Your VivaScan includes all software services, without you having to pay a subscription fee.

Intraoral scanning workflow



Consult

Show your patients their incredible new smile in just a few minutes with IvoSmile



Scan

Compact and intuitive intraoral scanning solution



Send to Lab

Fast and secure data transfer to your lab



Design & Produce in Lab

IPS e.max all-ceramic system for all indications of fixed prosthetics



Place

Permanent cementation of demanding ceramic and composite restorations



Care

Professional protection for your restoration