

Direct Restoratives



When **simplicity**
meets **esthetics**

Efficient
Esthetics


ivoclar
vivadent
passion vision innovation

Efficient Esthetics

Products that are optimally coordinated for direct restorative procedures enhance the efficiency with which high-quality esthetic restorations can be achieved.

Natural shade blend due to the chameleon effect of the Tetric Evo-Line composites

High bond strength (> 25 MPa) on wet and dry dentin thanks to Adhese® Universal

Reliable light-curing thanks to Bluephase® Style

Effective isolation and easier access due to OptraGate®

3 instead of 2 patients

treated in the same amount of time without compromises in quality¹

¹ Source: Special Feature "A fast and comfortable way of restoring teeth, Tetric EvoFlow® Bulk Fill / Tetric EvoCeram® Bulk Fill", Ivoclar Vivadent, January 2016

² A total of 890 fillings were monitored worldwide in seven clinical studies over observation periods ranging from 5 to 10 years. The overall survival rate calculated for Tetric EvoCeram® was 95.1 per cent.

³ Based on Tetric Evo-Line sales figures

A close-up photograph of two dental adhesive dispensers. The dispenser in the foreground is black with a green trigger and a clear nozzle. The background dispenser is also black with a green trigger. The background is a soft, out-of-focus white.

95%

survival rate²

More than

200

million restorations placed worldwide³

OptraGate®

Moisture control made easy

OptraGate® is a latex-free lip and cheek retractor. The aid's unique flexibility and elasticity make it comfortable to wear and assist patients to keep their mouth open.

Effective isolation
due to easier access

Focused and **efficient**
treatment procedures

Used in over

70

million treatments to date¹



Adhese[®] Universal

Universal adhesion with advanced delivery

Adhese[®] Universal is a single-component, light-cured adhesive for direct and indirect bonding procedures and all etching protocols. It is the only universal bonding agent available in the innovative VivaPen[®].

Universally suitable for direct and indirect bonding procedures and all etching protocols

Efficient VivaPen[®] delivery form provides almost 3 times more applications per millilitre compared to conventional bottles

Fast and convenient direct intraoral application with the ergonomic VivaPen[®]

High bond strength (> 25 MPa) on wet and dry dentin

Minimal risk of postoperative sensitivity due to the integrated desensitizing effect

**WET & DRY
BONDING**

Desensitizing
EFFECT

**Fast
and convenient**





190

single-tooth applications¹

Universal application
for direct and indirect
restorations

Tetric Evo-Line

One solution for all cavities

The Tetric Evo-Line range of products consists of four well-coordinated composite resins for your direct restorative needs.

Tetric EvoCeram® is a sculptable universal composite and **Tetric EvoFlow®** is a flowable universal composite for anterior and posterior restorations.

Tetric EvoCeram® Bulk Fill is a sculptable composite and **Tetric EvoFlow® Bulk Fill** a flowable composite for 4-mm bulk placement in the posterior region.

Advantages that convince users

- chameleon effect that ensures a natural shade blend
- true-to-nature esthetics for posterior restorations
- time savings of up to 47%¹
- one kit with only 6 materials for all cavities



More than

200

million restorations placed worldwide²

95%

survival rate³



¹ Tetric EvoCeram Bulk Fill and Tetric EvoFlow Bulk Fill compared to Tetric EvoCeram and Tetric EvoFlow; data available upon request.

² Based on Tetric Evo-Line sales figures

³ A total of 890 fillings were monitored worldwide in seven clinical studies over observation periods ranging from 5 to 10 years. The overall survival rate calculated for Tetric EvoCeram[®] was 95.1 per cent.

IPS Empress® Direct

The composite for highly esthetic fillings

IPS Empress® Direct is a highly esthetic composite for filling cavities in anterior and posterior teeth. It is available in five different translucency levels. Therefore, different types of anatomical features can be faithfully reproduced.

Lifelike shades for natural-looking restorations

Easy application of shade nuances with **Trans Opal**

Long-lasting esthetic results

Convenient handling for sophisticated requirements

Clever mastery of demanding restorations





More than

25

million restorations
placed worldwide¹

¹ Based on sales figures

Bluephase® Style

Reliable curing, short exposures

The innovative Bluephase® Style is designed for the intraoral polymerization of light-curing dental materials.

The Bluephase® curing lights are not only the best-selling LED polymerization lights in Europe¹, but they are also used in numerous clinical studies. They offer you the following advantages:

Reliable curing performance for
direct and indirect restorations

Short curing times starting from
5 seconds due to high light intensities²

Suitable for every light-curing material
due to polywave LED

Wide 10-mm light guide for time-saving
single-exposure curing procedures

1,200 mW/cm²



Find out more about the features that lend
uniqueness and high quality to every Bluephase:
www.ivoclarvivadent.com/bluephase_en





Simple two-button operation

Contactless battery charging

Ergonomically designed to comfortably fit every hand

Europe's number **1**¹
3-year warranty³

¹ Based on sales figures, acc. to KeyStone, 2017
² Only in conjunction with the Bluephase Style 20i
³ Warranty on battery: 1 year

Efficient esthetics – even better together

The combination of the Tetric Evo-Line with other state-of-the-art products from Ivoclar Vivadent allows the efficiency of the treatment procedure to be enhanced even further.

“I have analyzed the restorations placed in more than 100 patients based on before and after digital images with the date and time stamps. While I had to schedule 45 minutes for an MOD filling in the past, I only need 30 minutes today for this type of restoration in the same high quality and esthetics – thanks to Tetric EvoCeram® Bulk Fill, Adhese® Universal and the short polymerization times provided by Bluephase® Style.”

Dr Petr Hajný
Prague, Czech Republic




3 instead of 2

patients

treated in the same amount
of time without compromises
in quality

**Save valuable time during
appointments.**

The diagram compares the time scheduled
for conventional and efficiency-optimized
treatment protocols to place Class II
MOD restorations.



Source: Special Feature "A fast and comfortable way of restoring teeth, Tetric EvoFlow® Bulk Fill / Tetric EvoCeram® Bulk Fill", Ivoclar Vivadent, January 2016



More information:
www.ivoclarvivadent.com/hajny-en



1 Isolate

Effective isolation with
OptraGate®, cotton rolls and saliva
ejector



2 Bond

Direct intraoral application of
Adhese® Universal with the
efficient VivaPen®



3 Restore

Tetric Evo-Line – One solution
for all cavity classes

Lifelike, esthetic restorations
with **IPS Empress® Direct**



4 Contour

Time-saving contouring with
OptraSculpt® due to its anti-stick
effect



5 Cure

Reliable and fast curing with the
Bluephase® curing units



6 Polish

High-gloss polishing in only one step
with **OptraPol®**

Efficient Esthetics

Products that are optimally
coordinated for direct restora-
tive procedures enhance the
efficiency with which high-
quality esthetic restorations
can be achieved.

Production and distribution
Ivoclar Vivadent AG
Bendererstr. 2
9494 Schaan
Liechtenstein
Tel. +423 235 35 35
Fax +423 / 235 33 60
www.ivoclarvivadent.com

Distribution Germany
Ivoclar Vivadent GmbH
Dr. Adolf-Schneider-Str. 2
D-73479 Ellwangen, Jagst
Tel. +49 7961 8890
Fax +49 7961 6326
info@ivoclarvivadent.de
www.ivoclarvivadent.de

690976/DE/2018-09



ivoclar
vivadent[®]
passion vision innovation