

## OptraSculpt® Next Generation

OptraSculpt Next Generation is a contouring instrument for the efficient contouring of composite restorations in the posterior region. The two-components attachments produced in a special procedure facilitate the modelling of smooth and esthetic high-quality filling surfaces. The attachment shapes specifically designed for the posterior region allow for a natural-looking recreation of anatomical tooth structures such as cusp slopes, fissures and proximal margins.



### Indication:

Suitable for adapting, shaping and contouring uncured sculptable composite materials. Also suitable for lab composites.

### Advantages:


- Three attachment shapes for the professional creation of sophisticated anatomical tooth structures
- Optimum working angle for the posterior region
- Hygienic dispenser box for easy access to attachments

### Delivery forms

|   |             |               |
|---|-------------|---------------|
| <b>OptraSculpt® &amp; OptraSculpt® Pad System Kit</b> |             | <b>683066</b> |
| Instrument single-ended                               | 1 x         |               |
| Instrument double-ended                               | 1 x         |               |
| Pads 4 mm   | 60 x        |               |
| Ball  | 20 x        |               |
| Chisel  | 20 x        |               |
| Point   | 20 x        |               |
| <b>OptraSculpt® Starter Kit</b>                       |             | <b>683067</b> |
| Instrument double-ended                               | 1 x         |               |
| Ball  | 20 x        |               |
| Chisel  | 20 x        |               |
| Point   | 20 x        |               |
| <b>OptraSculpt® Refill ball</b>                       | <b>60 x</b> | <b>683068</b> |
| <b>OptraSculpt® Refill chisel</b>                     | <b>60 x</b> | <b>683069</b> |
| <b>OptraSculpt® Refill point</b>                      | <b>60 x</b> | <b>683100</b> |
| <b>OptraSculpt® Instrument double-ended Refill</b>    | <b>1 x</b>  | <b>683154</b> |

## OptraSculpt® Pad

OptraSculpt® Pad is a contouring instrument with special foam pad attachments designed for the efficient, non-stick forming and shaping of composites. The highly flexible synthetic foam pads optimally adjust to the anatomical contours and allow smooth modelling, without leaving any unwanted marks on the composite restoration.



Reference scale on the instrument handle:

The continuous vertical lines show ideal average widths, the dotted vertical lines average axial inclinations of the maxillary anterior teeth.

### Indication:

Suitable for adapting, shaping and contouring uncured sculptable composite materials. It is also available for lab composites.

### Advantages:

- Anti-stick effect: reduced adhesion of the composite to the instrument
- Professional esthetics due to homogeneous surface design and reference scales on the instrument handle
- Reduced effort required for finishing and polishing

### Delivery forms

|   |            |               |
|---|------------|---------------|
| <b>Assortment</b>                         |            | <b>639791</b> |
| OptraSculpt Pad Instrument                | 1x         |               |
| OptraSculpt Pad Refill/4 mm               | 60x        |               |
| OptraSculpt Pad Refill/6 mm               | 60x        |               |
| <b>OptraSculpt® Pad Refill/4 mm</b>       | <b>60x</b> | <b>639792</b> |
| <b>OptraSculpt® Pad Refill/6 mm</b>       | <b>60x</b> | <b>639793</b> |
| <b>OptraSculpt® Pad Instrument Refill</b> | <b>1x</b>  | <b>646004</b> |