

OptraGloss®

Efficient and advanced universal polishing system with a high concentration of diamond particles. A high-gloss finish and smooth surfaces are achieved in only one working step for composites and in only two steps for ceramics. OptraGloss polishers can be autoclaved and re-used.

Shapes:

- Flame PP
- Cup PP
- Lens PP
- Flame HP
- Cup HP
- Lens HP
- Spiral Wheel HP



Indications:

- Clinical polishing of glass-ceramic, feldspathic-ceramic and oxide-ceramic restorations
- Clinical polishing of composite resin restorations

Advantages:

- Radiant high-gloss restoration surfaces in only one step
- High wear resistance
- Autoclavable and reusable
- Two-step polishing system for ceramics

Delivery forms

Assortment		684548AN
Flame PP, Cup PP, Lens PP	2x each	
Flame HP, Cup HP, Lens HP, Spiral Wheel HP	2x each	
Ceramic Kit		684555AN
Flame PP, Cup PP, Lens PP	1x each	
Flame HP, Cup HP, Lens HP, Spiral Wheel HP	1x each	
Composite Kit		701870AN
Flame HP, Cup HP, Lens HP, Spiral Wheel HP	2x each	
Refill		
Flame PP	10x	684553AN
Cup PP	10x	684549AN
Lens PP	10x	684551AN
Flame HP	10x	684554AN
Cup HP	10x	684550AN
Lens HP	10x	684552AN
Spiral Wheel HP	10x	721094AN