

Astropol®

Indication:

Astropol is indicated for the pre-polishing and polishing of restorations made of composite material and amalgam.

Astropol® F

Polishing tips (grey) for the removal of excess composite material from proximal and occlusal surfaces

Astropol® P

Polishing tips (green) for polishing composite materials

Astropol® HP

Polishing tips (pink) for final high-gloss

Shapes:

- Small flame
- Large flame
- Cup
- Disc



Delivery forms

Assortment	557625
Astropol F: Small flame, large flame, cup, disc	3x each
Astropol P: Small flame, large flame, cup, disc	2x each
Astropol HP: Small flame, large flame, cup, disc	1x each

Refill Astropol® F (grey)

Small flame	6x	557615
Large flame	6x	557616
Cup	6x	557617
Disc	6x	557618

Refill Astropol® P (green)

Small flame	6x	557614
Large flame	6x	557621
Cup	6x	557619
Disc	6x	557620

Refill Astropol® HP (pink)

Small flame	3x	557613
Large flame	3x	557622
Cup	3x	557623
Disc	3x	557624

Politip®

Multiple-use, 2-step polishing system for composites and amalgam.

Indication:

Politip® is indicated for the polishing of restorations made of composite material and amalgam.

Shapes:

- Small cup
- Large cup
- Flame



Delivery forms

Assortment Politip® Polisher (P) / Finisher (F)	551580
Finisher shapes: flame	6x
large cup, small cup	3x each
Polisher shapes: flame	6x
large cup, small cup	3x each

Refill Politip® F

Flame	6x	533602
Large cup	6x	533603
Small cup	6x	533604

Refill Politip® P

Flame	6x	533584
Large cup	6x	533585
Small cup	6x	533587