

Tetric EvoFlow® Bulk Fill

Flowable and highly efficient 4-mm bulk-fill composite suitable as a first volume layer in Class I and II restorations.

Indications:

- As a first layer / first increment in composite restorations of Class I and II in permanent teeth
- Restoration of deciduous teeth

Tetric EvoFlow® Bulk Fill restorations are capped with e.g. Tetric EvoCeram® Bulk Fill.



Advantages:

Esthetics

- Chameleon effect ensures a natural shade blend
- Flowable dentin-like volume replacement material

Quality

- Layer thickness up to 4 mm
- Ideal combination with Tetric EvoCeram Bulk Fill

Efficiency

- Curing time of 10 seconds (> 1,000 mW/cm²)
- Time savings of up to 47% compared to the conventional technique*

*Compared with Tetric EvoCeram and Tetric EvoFlow; data available on request.

Shades:

Three dentin-like universal shades

- ^{IV}A: Universal A shade
- ^{IV}B: Universal B shade
- ^{IV}W: White shade for deciduous teeth or light-coloured permanent teeth

Delivery forms

Syringes

| Refill | 1 x 2 g | see p. 20 |
|--------------------------------|---------|-----------|
| 1 syringe, 2 g each (3 shades) | | |

Cavifils

| Refill | 20 x 0.2 g | see p. 20 |
|------------------------------------|------------|-----------|
| 20 Cavifils, 0.2 g each (3 shades) | | |

| | | |
|-----------------------------------|--------|--------|
| Shade guide | 1x | 667472 |
| Application tips ø 1.2 mm (black) | 1 x 20 | 669838 |

Tetric EvoCeram® Bulk Fill

Sculptable and highly efficient 4-mm bulk-fill composite for the posterior region.

Indications:

- Restoration of deciduous teeth
- Restoration of posterior teeth (Classes I and II, including reconstruction of individual cusps)
- Class V restorations (cervical caries, root erosion, wedge-shaped defects)
- Reconstructive build-up
- Extended fissure sealing in molars and premolars



Advantages:

Esthetics

- Chameleon effect ensures a natural shade blend
- Enamel-like volume replacement material

Quality

- Layer thickness up to 4 mm

Efficiency

- Curing time of 10 seconds (> 1,000 mW/cm²)
- Time savings of up to 47% compared to the conventional technique*

*Compared with Tetric EvoCeram and Tetric EvoFlow; data available on request.

Shades:

Three enamel-like universal shades

- ^{IV}A: Universelle A-Farbe
- ^{IV}B: Universelle B-Farbe
- ^{IV}W: White shade for deciduous teeth or light-coloured permanent teeth

Delivery forms

Syringes

| Refill | 1 x 3 g | see p. 20 |
|--|----------|-----------|
| 1 syringes, 3 g each (3 shades) | | |
| Refill | 3 x 3 g | see p. 20 |
| 3 syringes, 3 g each (3 shades) | | |
| Jumbo Refill | 10 x 3 g | 649792 |
| 10 syringes, 3 g each (^{IV} A) | | |

Cavifils

| Refill | 20 x 0.2 g | see p. 20 |
|---|-------------|-----------|
| 20 Cavifils, 0.2 g each (3 shades) | | |
| Jumbo Refill | 150 x 0.2 g | 649793 |
| 150 Cavifils, 0.2 g each (^{IV} A) | | |
| Shade guide | 1x | 667472 |