

## Tetric®

Light-curing, universal microhybrid composite with a special filler composition. Available in syringes and Cavifils.



before



after



before



after: with Tetric Basic White

### Indications:

- Anterior and posterior restorations

### Advantages:

- Excellent working properties
- Does not stick to the instrument
- Optimum modelling properties
- High abrasion resistance
- Very high radiopacity
- Fluoride release
- Natural shade effect (three degrees of translucency and ten shades available)
- Stability of shade

### Shades:

Universal shades: A1, A2, A3, A3.5, A4, B2, B3

Highly translucent enamel shade: T

Dentin shade: A3.5 Dentin

Basic White

### Delivery forms

#### Syringes

<b>Refill</b>	<b>1 x 3 g</b>	<b>see p. 21</b>
1 syringe, 3 g each (9 shades)		

<b>Refill</b>	<b>3 x 3 g</b>	<b>see p. 21</b>
3 syringes, 3 g each (Basic White)		

#### Cavifils

<b>Refill</b>	<b>20 x 0.2 g</b>	<b>see p. 21</b>
20 Cavifils, 0.2 g each (9 shades)		

<b>Refill</b>	<b>60 x 0.2 g</b>	<b>see p. 21</b>
60 Cavifils, 0.2 g each (Basic White)		

<b>Shade guide</b>	<b>1x</b>	<b>612524</b>
--------------------	-----------	---------------