

IPS Empress® Direct Color

For true-to-nature characterizations.



Indications:

- Individual characterization of direct and lab-fabricated composite restorations
- Masking of discoloured tooth structure
- Characterization of CAD/CAM veneers fabricated chairside

Advantages:

- True-to-nature characterization of high-end composite and CAD/CAM-fabricated restorations
- User-friendly direct application without dilution
- High masking capability even when applied in thin layers
- 0.4-mm tips for precise application
- The ideal complement to IPS Empress® Direct and IPS Empress CAD

Shades:

white, honey, ochre, brown, violet, blue, grey

Delivery forms

Assortment		637898
IPS Empress Direct Color	7 x 1 g each	
Shades: white, honey, ochre, brown, violet, blue, grey		
IPS Empress Direct Opaque	1 x 1.8 g	
var. accessories		

Refill	1 x 1 g	Possible applications	
white		<ul style="list-style-type: none"> • Imitating enamel stains • Simulating halo effects • Accentuating cusps and ridges 	637899
honey		<ul style="list-style-type: none"> • Emphasising the dentin core • Highlighting fissures • Imitating individual enamel characteristics 	637900
ochre			637903
blue		<ul style="list-style-type: none"> • Simulating translucent enamel areas • Imitating opalescent areas 	637901
violet			637902
grey			637905
brown		<ul style="list-style-type: none"> • Imitating discoloured fissures, pits and enamel cracks 	637904
Application tips ø 0.4 mm (black)			1 x 20 634847

IPS Empress® Direct Opaque

IPS Empress® Direct Opaque is a light-curing opaquing material suitable for masking exposed metal surface and discoloured tooth structures.



Indications:

- Masking of
- exposed metal surfaces when repairing damaged ceramic or composite veneers intraorally
 - discoloured tooth structure

Advantages:

- High masking capability even when applied in thin layers
- Universal opaque shade, no shade matching required

Shade:

opaque

Delivery forms

Assortment		637898
IPS Empress Direct Color	7 x 1 g each	
white, honey, ochre, brown, purple, blue, grey		
IPS Empress Direct Opaque	1 x 1.8 g	
var. accessories		
Refill	1 x 1.8 g	637897AN
1 syringe, 1.8 g (opaque)		
Application tips ø 0.9 mm (black)	1 x 20	634849