

Adhese® Universal

The universal adhesive

Adhese® Universal is a single-component, light-cured adhesive for direct and indirect bonding procedures that features compatibility with all etching techniques: self-etch, selective-enamel-etch and total-etch. Adhese Universal is available in the unique VivaPen® delivery for fast and convenient direct intra-oral application with minimal waste to optimize cost-effectiveness.



Indications:

- Direct-placed light-curing composite and compomer restorations
- Direct-placed core build-ups with light-, self- and dual-curing composites
- Adhesive cementation of indirect restorations with light- and dual-curing luting composites
- Repairs of fractured composite and compomer restorations
- Desensitization of hypersensitive cervical areas
- Sealing of prepared tooth surfaces before temporary/permanent cementation of indirect restorations

Advantages:

Esthetics

- Universally suitable for direct and indirect bonding procedures and all etching protocols

Quality

- High bond strength (> 25 MPa) on wet and dry dentin
- Minimal risk of postoperative sensitivity due to the integrated desensitizing effect

Efficiency

- Fast and convenient direct intraoral application with the ergonomic VivaPen
- Efficient VivaPen Delivery forms provides almost 3 times more applications per ml compared to conventional bottles*
- Curing time of only 3 seconds with the 3s PowerCure system: for restorations in the posterior region of permanent dentition (Class I and II) when light-cured from the occlusal aspect

* Data available on request

Delivery forms

VivaPen®

Refill	1 x 2 ml	665156
Refill	3 x 2 ml	664505
System Kit	1 x 2 ml	664504
Starter Kit	1 x 0.5 ml	663717
Refill VivaPen® snap-on cannulas	100x	627178
Refill VivaPen® snap-on cannulas	300x	627161
Refill VivaPen® protective sleeves	300x	627492

Bottle

Refill	1 x 5 g	663720
Refill	2 x 5 g	663721
Refill Vivadent applicators (green)	50x	664698
Refill Vivadent applicators (black)	50x	695128