

## ExcITE® F

ExcITE® F is a light-curing, fluoride-releasing single-component adhesive that is suitable for the total-etch technique.



### Indications:

- Direct light- and dual-cure composite and compomer restorations
- Adhesive luting of indirect restorations involving light-curing

### Advantages:

- Convenient and economical application with the VivaPen®
- Also available in bottles and single dose vessels
- High adhesion values on enamel and dentin
- Acetone-free (alcohol-based solvent)
- Fluoride release
- Reduced post-operative sensitivity
- Storage at room temperature

### Delivery forms

Intro Pack VivaPen®		627177
ExcITE F VivaPen	1 x 2 ml	
Total Etch Etching Gel	1 x 2 g	
var. accessories		
<b>Refill VivaPen, var. accessories</b>	<b>1 x 2 ml</b>	<b>613724</b>
<b>Refill VivaPen snap-on cannulas</b>	<b>100x</b>	<b>627178</b>
<b>Refill VivaPen snap-on cannulas</b>	<b>300x</b>	<b>627161</b>
<b>Refill VivaPen protective sleeves</b>	<b>300x</b>	<b>627492</b>
<b>Refill Bottle</b>	<b>1 x 5 g</b>	<b>630375</b>
<b>Refill Bottle</b>	<b>2 x 5 g</b>	<b>630376</b>
<b>Refill regular applicators (light blue)</b>	<b>50x</b>	<b>631699</b>
<b>Refill Single Dose</b>	<b>50 x 0.1 g</b>	<b>630377</b>

## ExcITE® F DSC

ExcITE® F DSC is a dual-cure, fluoride-releasing enamel and dentin adhesive that is suitable for the total-etch technique.



### Indications:

- Adhesive cementation of root posts
- Adhesive cementation of indirect, metal-free translucent restorations with dual- or self-curing luting composites
- Direct restorations and core build-ups with dual- or self-curing composites

### Advantages:

- Patented initiator-coated applicator brush
- Time-saving: one step, one coat
- No light-curing necessary
- High bonding values on enamel and dentin
- Single-dose presentation with applicators in two different sizes: Regular and Small / Endo
- Storage at room temperature

### Delivery forms

Refill		
Regular	50 x 0.1 g	630378
Small/Endo	50 x 0.1 g	630380