

Syntac®

Syntac® is the classic adhesive system for a chemically sound bond between the composite materials and the tooth structure.



Indications:

- Direct composite and compomer restorations in the anterior and posterior region
- Adhesive cementation of indirect restorations made of composite or ceramic

Advantages:

- For direct and indirect restorations
- Clinically proven for many years
- Reduction of post-operative sensitivity
- Excellent instant bond between dentin/enamel and composite
- High bond strength
- Storage at room temperature

Delivery forms

Assortment	2 x 3 g / 1x 6 g	532891
Syntac Primer, Syntac Adhesive, Heliobond, var. accessories		
Refill Primer	1 x 3 g	532893
Refill Adhesive	1 x 3 g	532892

Heliobond

Heliobond is a light-curing single-component bonding agent to optimize enamel etching with all light-curing restorative materials.



Indications:

- Adhesive for enamel bonding
- Together with Syntac also dentin bonding agent
- Cementation of tooth jewellery
- Transparent sealing of fissures and restorations
- As a bonding layer for the repair of composite, crown and bridge restorations

Advantages:

- For an optimum marginal seal between tooth and composite. Applied in conjunction with the enamel-etch technique
- Storage at room temperature

Delivery forms

Refill	1 x 11 g	532906
Refill	1 x 6 g	532907