

## Total Etch

Total Etch is a 37% phosphoric acid gel.



### Indications:

- Partial etching technique for adhesive restorations
- Total-etch technique for adhesive restorations
- Fissure sealings
- Cementation of brackets and tooth jewellery
- Adhesive cementation

### Advantages:

- The conventional etching gel with balanced handling properties
- Strong contrast to the tooth
- Pre-bent metal cannulas for accurate application
- Jumbo Delivery forms to keep adequate amounts of gel in stock

### Delivery forms

<b>Assortment Total Etch Jumbo</b>	<b>1 x 30 g / 2 x 2 g</b>	<b>558530</b>
Jumbo syringe application tips, empty 2-ml syringe		
<b>Refill</b>	<b>2 x 2 g</b>	<b>550588</b>
<b>Etch Refill application tips (blue)</b>	<b>1 x 20</b>	<b>550587</b>

## Email Preparator blue

37% phosphoric acid in a liquid Delivery forms.



### Indications:

- Enamel etching and dentin conditioning for sealants, composite restorations, and for the adhesive cementation of restorations.

### Advantages:

- Micro-interlocking of Heliobond bonding agent with the etched enamel
- The material flows well into the fissures due to the low viscosity
- Reliable retention of the sealant, composite, and the seated restoration

### Delivery forms

<b>Refill</b>	<b>1 x 6 g</b>	<b>533034</b>
---------------	----------------	---------------