

IPS® Natural Die Material Shade Guide

Determination of the preparation shade in the fabrication of all-ceramic restorations

The shade guide enables optimum selection of the shade after preparation. The shade is selected on the basis of the most severely discoloured part of the prepared tooth.

The IPS Natural Die Material Shade Guide allows you to fabricate a customized model die with which you can then accurately select the shade and brightness of all-ceramic restorations.



Indications:

Determination of the preparation shade after tooth preparation

Advantages:

- Easy and accurate reproduction of shades due to having determined the preparation shade
- Fabrication of a customized model die
- Autoclavability

Shades:

9 (ND 1, ND 2, ND 3, ND 4, ND 5, ND 6, ND 7, ND 8, ND 9)

Delivery forms

IPS® Natural Die Material Kit 597078

IPS Natural Die Material, 9 shades	1 x 8 g each
IPS Natural Die Material Separator	1 x 20 ml
IPS Natural Die Material shade guide	1x
var. accessories	

IPS® Natural Die Material Kit Promo 638177

IPS Natural Die Material Kit (597078)	1x
IPS Natural Die Material shade guide	4x

Refill IPS® Natural Die Material 1 x 8 g each

ND 1	597080
ND 2	597081
ND 3	597082
ND 4	597083
ND 5	597084
ND 6	597085
ND 7	597086
ND 8	597087
ND 9	597088

IPS® Natural Die Material Separator 1 x 20 ml **597089**

Die pins 10 x **554028**

IPS® Condenser 10 x **554030**

IPS® Natural Die Material shade guide (autoclavable) 1x **597090**

Example of the preparation shade effect



Restoration: Veneer (0.5 mm, IPS e.max® CAD HT B1)
Cementation material: Variolink® Esthetic neutral
Preparation shade: IPS® Natural Die Material, ND 1 – ND 9