

## Telio® CS C&B

Self-curing composite material in paste form for the fabrication of temporary restorations.



### Indications:

Fabrication of temporary restorations

- Crowns and bridges
- Inlays and onlays
- Core build-ups in combination with temporary posts
- Veneers

### Advantages:

- High fracture strength and optimum wear period
- Accuracy of fit and tension-free results
- Protection of the dental pulp and tissues
- Easy handling due to automix double-barrel cartridge

### Shades:

Bleach BL3, A1, A2, A3, A3.5

### Delivery forms

#### Telio® CS C&B Starter Kit A2 628244

Telio CS C&B A2	1 x 78 g
Mixing Tips	10 x
Telio CS Link Transparent	1 x 6 g
Mixing Tips	4 x
Telio CS Desensitizer	1 x 5 g
Application Brushes	50 x
var. accessories	

#### Refill Telio® CS C&B 1 x 78 g

Double-barrel cartridge, 10:1, var. accessories	
Bleach BL3	628249
A1	628245
A2	628246
A3	628247
A3.5	628248

#### Refill Dispenser 1x 561381

#### Refill Mixing Tips 50x 628157