

Telio® CS Link

Dual-curing, translucent, temporary composite cement for the esthetic and provisional cementation of temporary restorations.

Indications:

Temporary cementation of provisionals

- Crowns and bridges
- Inlays and onlays
- Veneers



Advantages:

- Eugenol-free
- Dual-curing
- High esthetics due to two translucent shades
- Ideal retention behaviour
- Easy handling due to automix double-push syringe

Shades:

A3 and Transparent

Delivery forms

| Refill Telio® CS Link | 2 x 6 g | |
|--|---------|--------|
| Double-push syringe, 1:1, var. accessories | | |
| A3 | | 627912 |
| Transparent | | 627914 |

| Refill Mixing Tips | 15x | 645951 |
|--------------------|-----|--------|
|--------------------|-----|--------|

| Telio® CS C&B Starter Kit A2 | | 628244 |
|------------------------------|--|--------|
|------------------------------|--|--------|

| | | |
|---------------------------|----------|--|
| Telio CS C&B A2 | 1 x 78 g | |
| Mixing Tips | 10 x | |
| Telio CS Link Transparent | 1 x 6 g | |
| Mixing Tips | 4 x | |
| Telio CS Desensitizer | 1 x 5 g | |
| Application Brushes | 50 x | |
| var. accessories | | |

Telio® CS Desensitizer

A clear solution for application onto dentin for reducing and preventing dentinal hypersensitivity and post-operative sensitivity.

Indications:

Prevention and reduction of dentinal hypersensitivity and post-operative sensitivity for the following fields of application:

- Temporary restorations
- Luting / cementation of indirect restorations
- Direct restorative treatment



Advantages:

- Total working time of 10 seconds. No light-curing required
- License for treatment planning and material selection. No interaction with the materials used in the further treatment steps.
- Reduction of dentinal hypersensitivity and post-operative sensitivities
- Cleaning of cut dentin surfaces

Delivery forms

| Refill Telio® CS Desensitizer | | |
|--|------------|--------|
| Refill Telio CS Desensitizer Bottle | 1 x 5 g | 627911 |
| Refill Telio CS Desensitizer Single Dose | 50 x 0.1 g | 627916 |

| Telio® CS C&B Starter Kit A2 | | 628244 |
|------------------------------|--|--------|
|------------------------------|--|--------|

| | | |
|---------------------------|----------|--|
| Telio CS C&B A2 | 1 x 78 g | |
| Mixing Tips | 10 x | |
| Telio CS Link Transparent | 1 x 6 g | |
| Mixing Tips | 4 x | |
| Telio CS Desensitizer | 1 x 5 g | |
| Application Brushes | 50 x | |
| var. accessories | | |