

Telio® CS Inlay

Telio® CS Inlay is a light-curing, single-component temporary restorative material.



Indications:

- Deep inlay preparations with parallel walls
- Reversible obturation of screw access holes in implant superstructures
- Relining of prefabricated temporary crowns made of polycarbonate or methacrylate

Advantages:

- Light-curing, single-component material
- Quick and easy application, removal in one piece
- No additional temporary cement or adhesive required
- Bonds to the tooth but does not stick to the instrument
- Demonstrates low polymerization shrinkage
- Seals tooth cavities effectively and shows no discolouration
- Available in two different shades and presentations (syringes and Cavifils)

Shades:

Transparent and Universal

Telio® CS Onlay

Telio® CS Onlay is a light-curing, single-component temporary restorative material.



Indications:

- Wide, shallow cavities with few retentive areas, such as onlay preparations

Advantages:

- Light-curing, single-component material
- Quick and easy application, removal in one piece
- No additional temporary cement or adhesive required
- Bonds to the tooth but does not stick to the instrument
- Demonstrates low polymerization shrinkage
- Seals tooth cavities effectively and shows no discolouration
- Available in two different shades and presentations (syringes and Cavifils)

Shades:

Transparent and Universal

Delivery forms

Refill Telio® CS Inlay Syringe	3 x 2.5 g	
Universal		635332
Transparent		635335

Refill Telio® CS Inlay Cavifil	30 x 0.25 g	
Universal		635333
Transparent		635334

Refill Telio® CS Onlay Syringe	3 x 2.5 g	
Universal		635339
Transparent		635338

Refill Telio® CS Onlay Cavifil	30 x 0.25 g	
Universal		635337
Transparent		635336