


SpeedCEM® Plus

SpeedCEM® Plus is a self-adhesive, self-curing resin cement with optional light-curing. The material offers an ideal combination of performance and user-friendliness. Its formulation has been optimized to make it particularly suitable for use in conjunction with restorations made of zirconium oxide and metal-ceramics and for the cementation of restorations on implant abutments.



SpeedCEM Plus is suited for the permanent seating of crowns and bridges on natural teeth and on implant abutments that are made of the following restorative materials:

- Oxide ceramics such as zirconium oxide, e.g. IPS e.max® ZirCAD
- Metal and metal-based restorations
- Lithium disilicate glass-ceramics, e.g. IPS e.max®

Advantages:

- Excellent self-cure performance
- User-friendly handling
- Easy clean-up
- Tolerant processing
- High radiopacity
- Efficient process

Shades:

Transparent, Yellow, Opaque

Delivery forms

Refills

SpeedCEM Plus Automix syringe	1 x 9 g	
Mixing tips	15x	
Root Canal Tips	5x	
Yellow		681710
Opaque		681711
Transparent		681712

Set of three Refills

SpeedCEM Plus Automix syringe	3 x 9 g	
Mixing tips	45x	
Root Canal Tips	15x	
Yellow		692423
Transparent		692422

Starter Pack

SpeedCEM Plus Automix syringe	1 x 2.5 g	
Ivoclean	1 g	
Mixing tips	4x	
Transparent		681713

System Kit

SpeedCEM Plus Automix syringe	1 x 9 g	
Ivoclean	5 g	
Mixing tips	15x	
Transparent		681714

Refill Mixing tips	15x	645951
Root Canal Tips	5x	645955