


Multilink® Automix

Multilink Automix® is an adhesive luting system for the cementation of indirect restorations made of silicate and oxide ceramics (e.g. IPS e.max®), metal, metal-ceramics as well as composites.



Indications:

Cementation of indirect restorations (inlays, onlays, crowns, bridges, endodontic posts) made of:

- Metal and metal-ceramics (PFM)
- All-ceramics (silicate)
- Reinforced all-ceramics (zirconium and aluminium oxide)
- Composites and fibre-reinforced composites (FRC)

Advantages:

- Strong hold – dual- as well as self-curing
- Universal – cementation of silicate and oxide ceramics as well as metal
- Clinically proven – many long-term studies
- Optimum removal of excess material
- Optimum esthetics
- Optimum marginal quality

Shades:

Transparent, yellow, white, opaque

Delivery forms

System Pack

| | | |
|-----------------------------------------|-----------|---------------|
| Multilink Automix syringe | 1 x 9 g | |
| Mixing Tips | 15 x | |
| Root Canal Tips | 5 x | |
| Multilink Primer A+B | 2 x 3 g | |
| Monobond Plus | 1 x 5 g | |
| Liquid Strip | 1 x 2.5 g | |
| Application Tips for Liquid Strip | 5 x | |
| Applicators (large and small) | 25 x each | |
| Mixing Pad | 1 x | |
| Mixing Plate (for Multilink Primer A+B) | 1 x | |
| var. accessories | | |
| Transparent | | 627471 |
| Yellow | | 627473 |
| Opaque | | 627472 |
| White | | 645954 |

Refill Multilink® Automix

| | | |
|---------------------------|---------|---------------|
| Multilink Automix syringe | 1 x 9 g | |
| Mixing Tips | 15x | |
| var. accessories | | |
| Transparent | | 615216 |
| Yellow | | 615217 |
| Opaque | | 615218 |
| White | | 645952 |

Refill Multilink® Automix Try-In

| | | |
|----------------------------------|-----------|---------------|
| Multilink Automix Try-In syringe | 1 x 1.7 g | |
| Transparent | | 645956 |
| Yellow | | 645957 |
| Opaque | | 645958 |
| White | | 645959 |

| | | |
|------------------------------------------------------|------------------|---------------|
| Refill Multilink® Primer A+B | 2 x 3 g | 576825 |
| Refill Multilink® Primer A | 1 x 3 g | 613626 |
| Refill Multilink® Primer B | 1 x 3 g | 613627 |
| Refill Monobond® Plus | 1 x 5 g | 626221 |
| Refill Liquid Strip | 1 x 2.5 g | 532505 |
| Refill Mixing Tips | 15x | 645951 |
| Root Canal Tips (Intra Canal and Mixing Tips) | 5 x each | 645955 |
| Application Tips for Liquid Strip | 20 x | 533666 |
| Refill Applicators (25 large, 25 small) | 2x 25 | 592435 |