

## Ivoclean®

Universal cleaning paste especially designed for removing contaminants, in particular those containing proteins, from restoration surfaces after intraoral try-in.



### Indication:

- Extraoral cleaning of pre-treated ceramic and metal restoration surfaces

### Delivery forms

Refill	1 x 5 g	637568
--------	---------	--------

## Liquid Strip

Glycerine gel to prevent the oxygen-inhibited layer of composites with composite or ceramic restorations.



### Indications:

- Try-in paste
- Preventing an inhibited layer with composites
- Protecting the etching pattern of lab-fabricated all-ceramic restorations

### Delivery forms

Refill	1 x 2.5 g	532505
--------	-----------	--------

## Heliosit® Orthodontic

Light-curing, highly transparent single-component luting material for brackets.



### Advantages:

- No dosing, no mixing
- Brackets can be placed on the tooth surfaces without time pressure due to the light-curing properties of the material

### Delivery forms

Refill	3 x 2.5 g	532629
--------	-----------	--------