

ApexCal®

ApexCal® is the new, creamy calcium hydroxide paste suitable for use as a temporary disinfectant dressing during endodontic treatment.



Direct application of ApexCal with a tip

Application of ApexCal with a lentulo spiral

Rinsing out of the temporary ApexCal dressing

Indications:

- Applied as a drug filler during root canal treatment and used for indirect and direct pulp capping

Advantages:

- Time-saving direct application
- Easy removal, as the material does not set when in place
- High radiopacity (400% Al)

Delivery forms

Refill		595915
Syringes ApexCal	2 x 2.5 g	
Cannulas	15x	
Refill Cannulas	15x	596608

Apexit® Plus

Calcium hydroxide root canal sealer for the permanent obturation of root canals.



Initial situation

Direct application of Apexit Plus using an Intra Canal Tip

Application of Apexit Plus using a lentulo spiral



Insertion of the gutta percha points

Condensation of the gutta percha points using a spreader

Obturated root canal

Indications:

- Permanent obturation of root canals in conjunction with all obturation techniques involving gutta-percha

Advantages:

- Excellent tissue tolerance
- Durable sealing of the root canal due to the slight setting expansion
- Its easy flowing composition allows the material to adapt well even to morphologically complicated canals
- Convenient application (automix syringe and Intra Canal Tip enable easy direct application)
- Long working time (mixed Apexit® Plus can be used over three hours at room temperature)

Delivery forms

Refill		593991
Double-push Syringes	2 x 6 g	
Mixing Tips	15x	
Intra Canal Tips	5x	
Refill Mixing Tips	15x	595132
Refill Intra Canal Tips	15x	597380