

FRC Postec® Plus

FRC Postec® Plus is a highly esthetic and highly radiopaque root canal post composed of glass fibre-reinforced composite. The post is suitable for the chairside reconstruction of endodontically treated teeth that show extensive coronal damage.

The FRC Postec Plus system consists of light-transmitting glass fibre-reinforced composite posts in three sizes and includes purpose-designed stainless steel reamers for the preparation of the recipient site of the post space.



Prepare the post space



Cement the post with Multilink



Build up the core with MultiCore Flow



Completed core build-up

Indications:

- Restoration of endodontically treated anterior and posterior teeth that demonstrate severe coronal damage

Advantages:

- Lifelike appearance
- Superior radiopacity
- High retention
- Gentle on roots
- Clearly arranged system using tried-and-tested products
- Easy to remove
- Clinically proven concept

Delivery forms

Intro Pack Size 0/1/3		603543
FRC Postec Plus Size 0	3x	
FRC Postec Plus Size 1	5x	
FRC Postec Plus Size 3	2x	
FRC Postec Plus Reamer Size 0	1x	
FRC Postec Plus Reamer Size 1	1x	
FRC Postec Plus Reamer Size 3	1x	
Refill Size 0	20x	603540
Refill Size 0	5x	603541
Refill Size 1	20x	590222
Refill Size 1	5x	583244
Refill Size 3	20x	590223
Refill Size 3	5x	583245
Refill Reamer, Size 0	1x	603542
Refill Reamer, Size 1	1x	572801
Refill Reamer, Size 3	1x	572800