

Virtual®

Virtual® is the new line of vinylpolysiloxane impression materials specifically developed to meet the requirements of today's most popular impression taking techniques.

The system features the five viscosities Putty, Heavy Body, Monophase, Light Body and Extra Light Body in both "fast" (2:30 min.) and "regular" (4:30 min.) setting times.

The system is complemented by a full line of accessories, including a tray adhesive.



Indications:

- Final impressions used for the fabrication of indirect restorations (crowns, bridges, inlays, onlays and veneers)
- Dental implant impressions
- Matrix for "wax-ups" and treatment planning, study models
- Edentulous impressions
- Matrix used to create temporary restorations

Advantages:

- Excellent adaptation to the moist dentin and soft tissues
- Excellent moisture displacement
- Exceptional readability and precision
- Flexible working times and fast setting phase
- Fresh peppermint scent
- Outstanding balance of physical properties

Shades:

- Sea blue: Virtual Putty, Virtual Heavy Body and Virtual Monophase
- Caramel: Virtual Light Body and Virtual Extra Light Body

Delivery forms

Multi Pack	12 x 50 ml + Mixing Tips	
Extra Light Body Fast		572812
Extra Light Body Regular		572813
Light Body Fast		572826
Light Body Regular		572827
Monophase Fast		572828
Monophase Regular		572829
Heavy Body Fast		572830
Heavy Body Regular		572831

Refill	2 x 50 ml + 12 Mixing Tips	
Extra Light Body Regular		562827
Extra Light Body Fast		562828
Light Body Regular		562829
Light Body Fast		562830

Refill	2 x 50 ml + 6 Mixing Tips	
Monophase Regular		562831
Monophase Fast		562832
Heavy Body Regular		562835
Heavy Body Fast		562836

Refill		
Putty Regular	2 x 300 ml	562839
Putty Fast	2 x 300 ml	562840

Accessories

Tray Adhesive	1 x 10ml	562841
Refill Mixing Tips, small	48x	572806
Refill Mixing Tips, large	48x	572807
Refill Intra Oral Tips, large	100x	572808
Refill Intra Oral Tips, small	100x	572809
Refill Dispenser 1:1	1x	572810