

## Virtual® 380

The two consistencies Heavy Body and Monophase are now offered as Virtual® 380 in new, large cartridges. These XXL cartridges are suitable for use in e.g. Pentamix, Plug & Press Dispenser and MixStar mixing units.



### Advantages:

- Ready-to-use, large cartridge - convenient handling
- Separate activation as in the case of foil bags is not required
- Economical package size - contents 380ml

### Delivery forms

<b>Multi Pack</b>		<b>6 x 380 ml each + Mixer + Bayonet Ring</b>
Monophase Fast		<b>594835</b>
Heavy Body Fast		<b>594836</b>
Heavy Body Regular		<b>594837</b>

<b>Refill</b>		<b>2 x 380 ml each + Mixer + Bayonet Ring</b>
Monophase Fast		<b>594838</b>
Heavy Body Fast		<b>594839</b>
Heavy Body Regular		<b>594840</b>

### Accessories

<b>Refill Dynamic Mixer</b>	<b>50x</b>	<b>600158</b>
-----------------------------	------------	---------------

Technical Data	Extra Light Body	Light Body	Monophase		Heavy Body		Putty
Colour	Caramel	Caramel	Sea blue		Sea blue		Sea blue
Contents	50 ml	50 ml	50 ml	380 ml	50 ml	380 ml	2 x 300 ml
Working time							
Fast	1:00 min.	1:00 min.	1:00 min.	1:00 min.	1:00 min.	1:00 min.	1:00 min.
Regular	1:30 min.	1:30 min.	1:30 min.	–	1:30 min.	1:30 min.	1:30 min.
Time in mouth							
Fast	2:30 min.	2:30 min.	2:30 min.	3:30 min.	2:30 min.	2:30 min.	2:30 min.
Regular	4:30 min.	4:30 min.	4:30 min.	–	4:30 min.	4:30 min.	4:30 min.
Linear dimensional change 24 h	≤ 1.50 %	≤ 1.50 %	≤ 1.50 %		≤ 1.50 %		≤ 1.50 %
Strain in compression (min.–max.)	2–8 %	2–8 %	2–8 %		0.8–5 %		0.8–5 %
Recovery from deformation	≥ 96.5 %	≥ 96.5 %	≥ 96.5 %		≥ 96.5 %		≥ 96.5 %
Detail reproduction	20 μ	20 μ	20 μ		50 μ		75 μ
Shore-A hardness	43–51	50–57	64–71		66–76		58–64

