

Silamat® S6

Form and function – perfectly mixed.

Indication:

- Modern universal mixer for amalgam, glass ionomer cements and other predosed dental materials in capsules and injection capsules



Advantages:

- Attractive design for the modern dental practice
- Thanks to its rounded shape and the removable cover the unit can be cleaned easily
- Tried-and-tested Silamat technology
- Above-average speed of 4500 rpm
- Special figure-eight pattern and patented mixing fork
- Compatible with all popular capsules and injection capsules in various shapes and sizes
- Three-year warranty
- Voltage 100–240 VAC

Delivery forms

Silamat® S6	100 – 240 VAC	602286
-------------	---------------	--------